

## F.I.R.S.T. in Iowa South Central [IASC] Soybean Entries - 2009

Company	Product*	Technology	SCN	Seed Trmt	MAT
<b>Channel</b>	3002R2*	RR2Y	R	AC	3.0
<b>Dyna-Gro</b>	33D27	RR	R	CM	2.7
	33RY30*	RR2Y	R	AC	3.0
	37A31	RR	R	CM	3.1
	37RY29*	RR2Y	R	AC	2.9
<b>FS Seeds</b>	A09-29	RR2Y	R	AC	2.9
	A09-301	RR2Y	R	AC	3.0
	A09-311	RR2Y	R	AC	3.1
	HS29R80	RR	R	CM	2.9
	HS30R72	RR	R	CM	3.0
<b>iCORN.com</b>	2.960R2*	RR2Y	R	AC	2.9
	3.070R2*	RR2Y	MR	AC	3.0
<b>Kruger</b>	K2-2601*	RR2Y	R	AC	2.6
	K2-2701*	RR2Y	R	AC	2.7
	K2-2801	RR2Y	R	AC	2.8
	K2-2901	RR2Y	R	AC	2.9
	K2-2902*	RR2Y	R	AC	2.9
	K2-3002	RR2Y	R	AC	3.0
	K2-3102	RR2Y	R	AC	3.1
	K2-3201	RR2Y	R	AC	3.2
<b>Pioneer</b>	93M11	RR	R	CM	3.1
	93Y01	RR	R	CM	3.0
<b>Prairie Brand</b>	PB-2439NRR2*	RR2Y	R	AC	2.4
	PB-2558NRR	RR	R	G,T	2.4
	PB-2828NRR2*	RR2Y	R	AC	2.8
	PB-2959NRR2*	RR2Y	R	AC	2.8
	PB-3039NRR2*	RR2Y	R	AC	3.0
<b>Producers</b>	2408R2	RR2Y	S	AC	2.4
	2808NR2	RR2Y	R	AC	2.8
	3008NR2	RR2Y	R	AC	3.0
<b>SOI</b>	2716NRR	RR	R	None	2.7
	2728NRR	RR	R	None	2.7
	2895RR	RR	S	None	2.8
	3136NRR	RR	R	None	3.1
<b>Stine</b>	27RA20	RR2Y	R	AC	2.7
	29RA20	RR2Y	R	AC	2.9
<b>Titan Pro</b>	29M3*	RR2Y	R	AC	2.9
	30M2*	RR2Y	R	AC	3.0
	31M3	RR2Y	R	None	3.1
	7272N	RR	R	None	2.7
<b>Trisler</b>	2504RR(CN)	RR	R	None	2.5
	29X0R2	RR2Y	R	AC	2.9
	31Z0R2	RR2Y	R	AC	3.1
	32X0R2	RR2Y	R	AC	3.2
	32Z0R2	RR2Y	R	AC	3.2

<b>Company</b>	<b>Product*</b>	<b>Technology</b>	<b>SCN</b>	<b>Seed Trmt</b>	<b>MAT</b>
----------------	-----------------	-------------------	------------	------------------	------------

\* **Product Identifiers:** ^ = discontinued product, \* = seed lot not commercially available at planting

**Technologies:** LL = Liberty Link, RR = Roundup Ready, RR2Y = Roundup Ready 2 Yield,

ST = Sulfonylurea Tolerant Soybean (STS), Lo Lin = Low Linolenic acid

**Soybean Cyst Nematode:** R = resistant, MR = moderate resistance, S = susceptible, -- = info not provided

**Seed Treatments:** A - Allegiance, AM - ApronMaxx, AP - ApronXL, CM - CruiserMaxx, E - Excaliber, G - Gaucho,

PO - PunchOut, SG - Sure Gro, T - Trilex, T6 - Trilex 6000, -- = info not provided